


## FADS1 Antibody / Fatty acid desaturase 1 / Delta-5 desaturase (R35401)

Catalog No.	Formulation	Size
R35401-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3992
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This FADS1 antibody is available for research use only.



Western blot testing of human lung lysate with FADS1 antibody at 0.3ug/ml. Predicted molecular weight: ~52 kDa.

### Description

Additional name(s) for this target protein: Fatty acid desaturase 1

### Application Notes

Optimal dilution of the FADS1 antibody should be determined by the researcher.

### Immunogen

Amino acids QPSFEPTKNKELTDE were used as the immunogen for this FADS1 antibody.

## **Storage**

Aliquot and store the FADS1 antibody at -20oC.