


## FANCE Antibody (R33797)

Catalog No.	Formulation	Size
R33797-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2178
Applications	Immunohistochemistry (FFPE) : 3-6ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This FANCE antibody is available for research use only.



IHC testing of FFPE human placenta with FANCE antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

## Description

Additional name(s) for this target protein: Fanconi anemia, complementation group E

## Application Notes

Optimal dilution of the FANCE antibody should be determined by the researcher.

## Immunogen

Amino acids EHKSLADGGS were used as the immunogen for this FANCE antibody.

## **Storage**

Aliquot and store the FANCE antibody at -20oC.