


FANCE Antibody (R33797)

Catalog No.	Formulation	Size
R33797-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2178
Applications	Immunohistochemistry (FFPE) : 3-6ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This FANCE antibody is available for research use only.



IHC testing of FFPE human placenta with FANCE antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Fanconi anemia, complementation group E

Application Notes

Optimal dilution of the FANCE antibody should be determined by the researcher.

Immunogen

Amino acids EHKSLADGGS were used as the immunogen for this FANCE antibody.

Storage

Aliquot and store the FANCE antibody at -20oC.