


## FANCF Antibody (R35716)

Catalog No.	Formulation	Size
R35716-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	2188
<b>Applications</b>	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
<b>Limitations</b>	This FANCF antibody is available for research use only.



Western blot of HEK293 lysate overexpressing human FANCF-MYC probed with FANCF antibody [1ug/ml] in Lane A and anti-MYC [1/1000] in lane C. Mock-transfected HEK293 probed with FANCF antibody [1ug/ml] in Lane B. Predicted molecular weight: ~42kDa.

### Description

Additional name(s) for this target protein: Fanconi anemia, complementation group F

### Application Notes

Optimal dilution of the FANCF antibody should be determined by the researcher.

### Immunogen

Amino acids RPQEELEPGIHK were used as the immunogen for this FANCF antibody.

## **Storage**

Aliquot and store the FANCF antibody at -20oC.