


FANCL Antibody (R35379)

Catalog No.	Formulation	Size
R35379-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	55120
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This FANCL antibody is available for research use only.



Western blot testing of human bone marrow lysate with FANCL antibody at 1ug/ml.
Predicted molecular weight: ~43kDa.

Description

Additional name(s) for this target protein: Fanconi anemia, complementation group L; FANCL

Application Notes

Optimal dilution of the FANCL antibody should be determined by the researcher.

Immunogen

Amino acids QNLKDVLEIDFP were used as the immunogen for this FANCL antibody.

Storage

Aliquot and store the FANCL antibody at -20oC.