


Fas Antibody / CD95 (R34661)

Catalog No.	Formulation	Size
R34661-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	355
Localization	Cytoplasmic, membranous
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Fas antibody is available for research use only.



Western blot testing of MOLT4 lysate with Fas antibody at 0.1ug/ml. Predicted molecular weight: ~38 kDa (unmodified), 40-55 kDa (glycosylated).

Description

Additional name(s) for this target protein: TNF receptor superfamily, member 6, TNFRSF6, CD95

Application Notes

Optimal dilution of the Fas antibody should be determined by the researcher.

1. This FAS antibody will detect isoforms 1, 2 and 3.

Immunogen

Amino acids KTCRKHRKENQGSH were used as the immunogen for this Fas antibody.

Storage

Aliquot and store the Fas antibody at -20oC.