


Fascin Antibody (R35572)

Catalog No.	Formulation	Size
R35572-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6624
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This Fascin antibody is available for research use only.



Western blot testing of human skeletal muscle lysate with Fascin antibody at 1ug/ml.
Predicted molecular weight ~55kDa.

Description

Additional name(s) for this target protein: Fascin-1, FSCN1

Application Notes

Optimal dilution of the Fascin antibody should be determined by the researcher.

Immunogen

Amino acids RDTKKCAFRTHT were used as the immunogen for this Fascin antibody.

Storage

Aliquot and store the Fascin antibody at -20oC.