


FBXO4 Antibody (R35165)

Catalog No.	Formulation	Size
R35165-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	26272
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This FBXO4 antibody is available for research use only.



Western blot testing of human testis lysate with FBXO4 antibody at 0.5ug/ml. Predicted molecular weight: ~35kDa.

Description

Additional name(s) for this target protein: F-box protein 4

Application Notes

Optimal dilution of the FBXO4 antibody should be determined by the researcher.

1. This FBXO4 antibody is specific for isoforms 1 and 2.

Immunogen

Amino acids PTAGLPQRQIDG were used as the immunogen for this FBXO4 antibody.

Storage

Aliquot and store the FBXO4 antibody at -20oC.