

FBXO43 Antibody (R35015)

Catalog No.	Formulation	Size
R35015-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	286151
Applications	Western Blot : 1-3ug/ml Flow Cytometry : suitable Immunofluorescence : 5-10ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This FBXO43 antibody is available for research use only.



Description

Additional name(s) for this target protein: F-box protein 43, early mitotic inhibitor 2, EMI2

Application Notes

Optimal dilution of the FBXO43 antibody should be determined by the researcher.

Immunogen

Amino acids PGSAQSKRNLKRL were used as the immunogen for this FBXO43 antibody.

Storage

Aliquot and store the FBXO43 antibody at -20oC.