


FE65 Antibody (R34133)

Catalog No.	Formulation	Size
R34133-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	322
Localization	Cytoplasmic, membrane
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This FE65 antibody is available for research use only.



Western blot testing of NIH3T3 lysate with FE65 antibody at 0.1ug/ml. Predicted molecular weight: ~77kDa; the 97kDa isoform can be cleaved to yield a 65kDa C-terminal form.

Description

Additional name(s) for this target protein: Amyloid beta (A4) precursor protein-binding, family B, member 1; APBB1

Application Notes

Optimal dilution of the FE65 antibody should be determined by the researcher.

Immunogen

Amino acids GSLKPKRLGAHTP were used as the immunogen for this FE65 antibody.

Storage

Aliquot and store the FE65 antibody at -20oC.

References (2)