


Ferd3l Antibody (R35263)

Catalog No.	Formulation	Size
R35263-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	114712
Applications	Immunohistochemistry (FFPE) : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Ferd3l antibody is available for research use only.



IHC testing of a PFA-perfused cryosection of mouse embryo E13.5 with Ferd3l antibody at 1ug/ml (Alexa Fluor 488 detection, green). Nuclear counterstain with DAPI (blue).

Description

Additional name(s) for this target protein: Fer3-like

Application Notes

Optimal dilution of the Ferd3l antibody should be determined by the researcher.

Immunogen

Amino acids QFDERYQEVEGDE were used as the immunogen for this Ferd3l antibody.

Storage

Aliquot and store the Ferd3l antibody at -20oC.