


## FGF13 Antibody (R34447)

Catalog No.	Formulation	Size
R34447-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2258
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This FGF13 antibody is available for research use only.



Western blot testing of mouse brain lysate with FGF13 antibody at 0.5ug/ml. Predicted molecular weight: ~28kDa.

### Description

Additional name(s) for this target protein: Fibroblast growth factor 13

### Application Notes

Optimal dilution of the FGF13 antibody should be determined by the researcher.

1. This FGF13 antibody is specific for isoform 1.

## **Immunogen**

Amino acids KGKTSCDKNKLVFS were used as the immunogen for this FGF13 antibody.

## **Storage**

Aliquot and store the FGF13 antibody at -20oC.