


Fgf14 Antibody (R35599)

Catalog No.	Formulation	Size
R35599-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2259
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Fgf14 antibody is available for research use only.



Western blot testing of mouse brain lysate with Fgf14 antibody at 1ug/ml. Predicted molecular weight: ~28kDa.

Description

Additional name(s) for this target protein: Fibroblast growth factor 14

Application Notes

Optimal dilution of the Fgf14 antibody should be determined by the researcher.

1. This Syndecan 1 antibody will detect isoforms human 1a and mouse a.

Immunogen

Amino acids REQHWDRPSASR were used as the immunogen for this Fgf14 antibody.

Storage

Aliquot and store the Fgf14 antibody at -20oC.