


FGF23 Antibody (R35632)

Catalog No.	Formulation	Size
R35632-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8074
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This FGF23 antibody is available for research use only.



Western blot testing of human hippocampus lysate with FGF23 antibody at 0.3ug/ml.
Predicted molecular weight 28~32kDa.

Description

Additional name(s) for this target protein: Fibroblast growth factor 23

Application Notes

Optimal dilution of the FGF23 antibody should be determined by the researcher.

Immunogen

Amino acids RHTRSAEDDSERD were used as the immunogen for this FGF23 antibody.

Storage

Aliquot and store the FGF23 antibody at -20oC.