


FGF5 Antibody (R35068)

Catalog No.	Formulation	Size
R35068-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2250
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This FGF5 antibody is available for research use only.



Western blot testing of human placenta lysate with FGF5 antibody at 0.1ug/ml. Predicted molecular weight: ~30kDa.

Description

Additional name(s) for this target protein: Fibroblast growth factor 5

Application Notes

Optimal dilution of the FGF5 antibody should be determined by the researcher.

1. This FGF5 antibody is specific for isoform 1.

Immunogen

Amino acids PLSAPRKNTNSVK were used as the immunogen for this FGF5 antibody.

Storage

Aliquot and store the FGF5 antibody at -20oC.