


## FGF5 Antibody (R35068)

Catalog No.	Formulation	Size
R35068-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2250
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This FGF5 antibody is available for research use only.



Western blot testing of human placenta lysate with FGF5 antibody at 0.1ug/ml. Predicted molecular weight: ~30kDa.

## Description

Additional name(s) for this target protein: Fibroblast growth factor 5

## Application Notes

Optimal dilution of the FGF5 antibody should be determined by the researcher.

1. This FGF5 antibody is specific for isoform 1.

## Immunogen

Amino acids PLSAPRKNTNSVK were used as the immunogen for this FGF5 antibody.

## Storage

Aliquot and store the FGF5 antibody at -20oC.