


FGF8 Antibody (R34310)

Catalog No.	Formulation	Size
R34310-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2253
Localization	Secreted
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This FGF8 antibody is available for research use only.



Description

Additional name(s) for this target protein: Fibroblast growth factor 8

Application Notes

Optimal dilution of the FGF8 antibody should be determined by the researcher.

Immunogen

Amino acids HTTEQSLRFEFLNYP were used as the immunogen for this FGF8 antibody.

Storage

Aliquot and store the FGF8 antibody at -20oC.