


FGFR1 Antibody (R35549)

Catalog No.	Formulation	Size
R35549-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2260
Applications	Western Blot : 1-3ug/ml (1) ELISA (peptide) LOD : 1:8000
Limitations	This FGFR1 antibody is available for research use only.



Western blot testing of A549 lysate with FGFR1 antibody at 1ug/ml. Predicted molecular weight of detectable isoforms: 89~93kDa.

Description

Additional name(s) for this target protein: Fibroblast growth factor receptor 1, FGFR

Application Notes

Optimal dilution of the FGFR1 antibody should be determined by the researcher.

1. This Ab is specific for isoforms 1, 2, 10, 11, 14, ranging in size from 819 AA to 853 AA.

Immunogen

Amino acids RCRLRDDVQSIN were used as the immunogen for this FGFR1 antibody.

Storage

Aliquot and store the FGFR1 antibody at -20oC.