


FGL1 Antibody (R34028)

Catalog No.	Formulation	Size
R34028-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2267
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This FGL1 antibody is available for research use only.



Western blot testing of human liver lysate with FGL1 antibody at 0.1ug/ml. Predicted molecular weight: ~36kDa.

Description

Additional name(s) for this target protein: Fibrinogen-like 1, Hepassocin

Application Notes

Optimal dilution of the FGL1 antibody should be determined by the researcher.

Immunogen

Amino acids ETRVKQQQVKIKQ were used as the immunogen for this FGL1 antibody.

Storage

Aliquot and store the FGL1 antibody at -20oC.