


## FGR Antibody (R36277)

Catalog No.	Formulation	Size
R36277-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Dog, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2268
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This FGR antibody is available for research use only.



Western blot testing of mouse spleen lysate with FGR antibody at 0.5ug/ml. Predicted molecular weight: ~59kDa, routinely observed at 55~59kDa.

## Description

Additional name(s) for this target protein: Gardner Rasheed feline sarcoma viral (v fgr); Proto oncogene tyrosine protein kinase FGR

## Application Notes

Optimal dilution of the FGR antibody should be determined by the researcher.

## **Immunogen**

Amino acids TSAEPQYQPGDQT were used as the immunogen for this FGR antibody.

## **Storage**

Aliquot and store the FGR antibody at -20oC.