


## FGR Antibody (R36277)

Catalog No.	Formulation	Size
R36277-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Dog, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2268
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This FGR antibody is available for research use only.



Western blot testing of mouse spleen lysate with FGR antibody at 0.5ug/ml. Predicted molecular weight: ~59kDa, routinely observed at 55~59kDa.

## Description

Additional name(s) for this target protein: Gardner Rasheed feline sarcoma viral (v fgr); Proto oncogene tyrosine protein kinase FGR

## Application Notes

Optimal dilution of the FGR antibody should be determined by the researcher.

## Immunogen

Amino acids TSAEPQYQPGDQT were used as the immunogen for this FGR antibody.

## Storage

Aliquot and store the FGR antibody at -20oC.