


## Fibrillarin Antibody (R34637)

Catalog No.	Formulation	Size
R34637-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2091
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Fibrillarin antibody is available for research use only.



Western blot testing of HEK293 lysate with Fibrillarin antibody at 1ug/ml. Predicted molecular weight ~34kDa.

## Description

Additional name(s) for this target protein: FIB, FLRN; FBL

## Application Notes

Optimal dilution of the Fibrillarin antibody should be determined by the researcher.

## Immunogen

Amino acids KRVSISEGDDKIEYR were used as the immunogen for this Fibrillarin antibody.

## Storage

Aliquot and store the Fibrillarin antibody at -20oC.