


Fibrinogen alpha chain Antibody (R36088)

Catalog No.	Formulation	Size
R36088-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2243
Applications	Western Blot : 0.01-0.03ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Fibrinogen alpha chain antibody is available for research use only.



Western blot testing of platelet lysate with Fibrinogen alpha chain antibody at 0.01ug/ml. Predicted molecular weight ~70 kDa (alpha isoform), ~95 kDa (alpha-E isoform).

Description

Additional name(s) for this target protein: FGA, Fib2

Application Notes

Optimal dilution of the Fibrinogen alpha chain antibody should be determined by the researcher.

Immunogen

Amino acids STSYNRGDSTFES were used as the immunogen for this Fibrinogen alpha chain antibody. This sequence is common to the alpha and alpha-E isoforms.

Storage

Aliquot and store the Fibrinogen alpha chain antibody at -20oC.