


Filamin A Antibody (R33486)

Catalog No.	Formulation	Size
R33486-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2316
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Filamin A antibody is available for research use only.



Western blot testing of human uterus lysate with Filamin A antibody at 0.1ug/ml.
Predicted molecular weight ~281kDa.

Description

Additional name(s) for this target protein: FLNA

Application Notes

Optimal dilution of the Filamin A antibody should be determined by the researcher.

Immunogen

Amino acids DIRDAPQDFHPDRVK were used as the immunogen for this Filamin A antibody.

Storage

Aliquot and store the Filamin A antibody at -20oC.