


FIZZ1 Antibody (R33649)

Catalog No.	Formulation	Size
R33649-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	84666
Applications	Immunohistochemistry (FFPE) : 2-4ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This FIZZ1 antibody is available for research use only.



IHC testing of FFPE human colon with FIZZ1 antibody at 2.5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Resistin like beta; RETNLB

Application Notes

Optimal dilution of the FIZZ1 antibody should be determined by the researcher.

Immunogen

Amino acids DSVMDKKIKDVLNS were used as the immunogen for this FIZZ1 antibody.

Storage

Aliquot and store the FIZZ1 antibody at -20oC.