


FLAP Antibody (R34221)


Catalog No.	Formulation	Size
R34221-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	241
Applications	Western Blot : 0.1-0.3ug/ml Immunohistochemistry (FFPE) : 3-6ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This FLAP antibody is available for research use only.



Western blot testing of human U937 cell lysate with FLAP antibody at 0.1ug/ml.
Predicted molecular weight: ~18kDa.



IHC staining of FFPE human spleen with FLAP antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Five-lipoxygenase activating protein; ALOX5AP

Application Notes

Optimal dilution of the FLAP antibody should be determined by the researcher.

Immunogen

Amino acids HKVEHESRTQNGR were used as the immunogen for this FLAP antibody.

Storage

Aliquot and store the FLAP antibody at -20oC.