



# FLAP Antibody (R34221)

Catalog No.	Formulation	Size
R34221-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	241
Applications	Western Blot : 0.1-0.3ug/ml Immunohistochemistry (FFPE) : 3-6ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This FLAP antibody is available for research use only.





IHC staining of FFPE human spleen with FLAP antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

## **Description**

Additional name(s) for this target protein: Five-lipoxygenase activating protein; ALOX5AP

## **Application Notes**

Optimal dilution of the FLAP antibody should be determined by the researcher.

## **Immunogen**

Amino acids HKVEHESRTQNGR were used as the immunogen for this FLAP antibody.

## **Storage**

Aliquot and store the FLAP antibody at -20oC.