


## Flotillin 1 Antibody (Biotin Conjugate) (R35932BTN)

Catalog No.	Formulation	Size
R35932BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Pig
Format	Biotin Conjugate
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10211
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This Flotillin 1 antibody is available for research use only.



Western blot testing of K562 lysate with biotinylated Flotillin 1 antibody at 1ug/ml.  
Expected/observed molecular weight ~49 kDa.

### Description

Additional name(s) for this target protein: FLOT1

### Application Notes

Optimal dilution of the Flotillin 1 antibody should be determined by the researcher.

### Immunogen

Amino acids SISQVNHKPLRTA were used as the immunogen for this Flotillin 1 antibody.

## **Storage**

Aliquot and store the Flotillin 1 antibody at -20oC.