

FOLR1 Antibody / FR alpha [clone FOLR1/13424R] (V5900)

Catalog No.	Formulation	Size
V5900-100UG	0.2 mg/ml in 1X PBS with 0.05% BSA, 0.05% sodium azide	100 ug
V5900-20UG	0.2 mg/ml in 1X PBS with 0.05% BSA, 0.05% sodium azide	20 ug
V5900SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Species Reactivity	Guinea Pig, Human, Mouse, Rat
Format	Purified
Host	Rabbit
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG, kappa
Clone Name	FOLR1/13424R
UniProt	P15328
Localization	Cell membrane, Cytoplasm, Secreted
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml Western Blot : 2-4ug/ml
Limitations	This FOLR1/FR alpha antibody is available for research use only.

Application Notes

Optimal dilution of the FOLR1/FR alpha antibody should be determined by the researcher.

Immunogen

A recombinant fragment (around amino acids 1-200) of human FOLR1 protein corresponding to the extracellular domain (exact sequence is proprietary) was used as the immunogen for the FOLR1/FR alpha antibody.

Storage

FOLR1/FR alpha antibody with sodium azide - store at 2 to 8°C; antibody without sodium azide - store at -20 to -80°C.

