


## FOXA1 Antibody (F44167)

Catalog No.	Formulation	Size
F44167-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F44167-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Antigen affinity purified
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit Ig
<b>Purity</b>	Antigen affinity
<b>UniProt</b>	P55317
<b>Applications</b>	Western Blot : 1:1000
<b>Limitations</b>	This FOXA1 antibody is available for research use only.



## Description

This gene encodes a member of the forkhead class of DNA-binding proteins. These hepatocyte nuclear factors are transcriptional activators for liver-specific transcripts such as albumin and transthyretin, and they also interact with chromatin. Similar family members in mice have roles in the regulation of metabolism and in the differentiation of the pancreas and liver.

## Application Notes

Titration of the FOXA1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## **Immunogen**

A portion of amino acids 270-298 from the human protein was used as the immunogen for this FOXA1 antibody.

## **Storage**

Aliquot the FOXA1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.