


## FOXA1 Antibody [clone FOXA1/1519] (V3370)


Catalog No.	Formulation	Size
V3370-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	100 ug
V3370-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	20 ug
V3370SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug

[Bulk quote request](#)


<b>Species Reactivity</b>	Human, Rat
<b>Format</b>	Purified
<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal (mouse origin)
<b>Isotype</b>	Mouse IgG2a, kappa
<b>Clone Name</b>	FOXA1/1519
<b>Purity</b>	Protein G affinity chromatography
<b>Buffer</b>	1X PBS, pH 7.4
<b>UniProt</b>	P55317
<b>Gene ID</b>	3169
<b>Localization</b>	Nuclear
<b>Applications</b>	Western Blot : 1-2ug/ml Immunohistochemistry (FFPE) : 1-2ug/ml for 30 min at RT Immunofluorescence : 2-4ug/ml Flow Cytometry : 1-2ug/million cells
<b>Limitations</b>	This FOXA1 antibody is available for research use only.




IHC testing of FFPE human prostate carcinoma with FOXA1 antibody (clone FOXA1/1519). HIER: boil sections in 10mM Tris with 1mM EDTA, pH9, for 10-20 min.




IHC testing of FFPE rat colon with FOXA1 antibody (clone FOXA1/1519). HIER: boil sections in 10mM Tris with 1mM EDTA, pH9, for 10-20 min.




Immunofluorescent staining of PFA-fixed human MCF7 cells with FOXA1 antibody (clone FOXA1/1519, green) and Phalloiden (red).



Western blot testing of human HeLa cell lysate with FOXA1 antibody (clone FOXA1/1519). Predicted molecular weight: ~49 kDa.



FACS testing of PFA-fixed human MCF-7 cells with FOXA1 antibody (clone FOXA1/1519); Red=isotype control, Blue= FOXA1 antibody.



SDS-PAGE Analysis of Purified, BSA-Free FOXA1 Antibody (clone FOXA1/1519). Confirmation of Integrity and Purity of the Antibody.

## Description

The transcription factor Forkhead-box A1 (FOXA1), also known as hepatocyte nuclear factor 3-alpha, is a member of the FOX class of transcription factors. HNF-1 (HNF-1 $\alpha$  and HNF-1 $\beta$ ), HNF-3 (HNF-3 $\alpha$ , HNF-3 $\beta$  and HNF-3 $\gamma$ ), HNF-4 (HNF-4 $\alpha$  and HNF-4 $\beta$ ), and HNF-6 compose, in part, a homeoprotein family designated the hepatocyte nuclear factor family. The various HNF-1 isoforms regulate transcription of genes in the liver as well as in other tissues such as kidney, small intestine and thymus. FOXA1 is expressed in normal breast ductal epithelium and other epithelium in different organs, such as lung, pancreas, bladder, prostate, and colon. Recently, FOXA1 has been shown to be a major determinant of estrogen-ER activity and endocrine response in breast cancer cells. FOXA1 expression correlates with estrogen receptor (ER)-positivity, especially

in luminal subtype A breast cancers, which is associated with favorable prognosis. FOXA1 is useful in the sub-classification of breast carcinomas.

## Application Notes

The concentration stated for each application is a general starting point. Variations in protocols, secondaries and substrates may require the FOXA1 antibody to be titered up or down for optimal performance.

## Immunogen

A human partial recombinant protein corresponding to amino acids 372-472 was used as the immunogen for this FOXA1 antibody.

## Storage

Store the FOXA1 antibody at 2-8oC (with azide) or aliquot and store at -20oC or colder (without azide).

## References (1)