


FOXB1 Antibody (R36293)

Catalog No.	Formulation	Size
R36293-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Pig, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	27023
Applications	Western Blot : 1-3ug/ml IF/ICC : suitable ELISA (peptide) LOD : 1:32000
Limitations	This FOXB1 antibody is available for research use only.



Western blot testing of HepG2 nuclear lysate with FOXB1 antibody at 1ug/ml. Predicted molecular weight: ~35kDa.

Description

Additional name(s) for this target protein: Forkhead box B1; FKH5

Application Notes

Optimal dilution of the FOXB1 antibody should be determined by the researcher.

Immunogen

Amino acids TSPASALHSVAVH were used as the immunogen for this FOXB1 antibody.

Storage

Aliquot and store the FOXB1 antibody at -20oC.