


FOXC1 Antibody (R33387)

Catalog No.	Formulation	Size
R33387-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2296
Applications	Western Blot : 2-4ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This FOXC1 antibody is available for research use only.



Western blot testing of HEK293 lysate with FOXC1 antibody at 2ug/ml. Predicted molecular weight: ~57 kDa, routinely observed at up to 75 kDa.

Description

Additional name(s) for this target protein: Forkhead box C1

Application Notes

Optimal dilution of the FOXC1 antibody should be determined by the researcher.

Immunogen

Amino acids DAVKDKEEKDRLH were used as the immunogen for this FOXC1 antibody.

Storage

Aliquot and store the FOXC1 antibody at -20oC.