


Foxi3 Antibody (R35677)

Catalog No.	Formulation	Size
R35677-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	232077
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This Foxi3 antibody is available for research use only.



Western blot testing of mouse skin lysate with Foxi3 antibody at 1ug/ml. Predicted molecular weight: ~42kDa.

Description

Additional name(s) for this target protein: Forkhead box i3

Application Notes

Optimal dilution of the Foxi3 antibody should be determined by the researcher.

Immunogen

Amino acids RLRVSGKLEGDSPA were used as the immunogen for this Foxi3 antibody.

Storage

Aliquot and store the Foxi3 antibody at -20°C.