


FOXK2 Antibody (R36241)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R36241-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|---------------------------|--|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Format | Antigen affinity purified |
| Host | Goat |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 3607 |
| Applications | Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000 |
| Limitations | This FOXK2 antibody is available for research use only. |



Western blot testing of human thymus lysate with FOXK2 antibody at 0.3ug/ml.
Predicted molecular weight: ~69kDa (isoform 1).

Description

Additional name(s) for this target protein: Forkhead box K2; interleukin enhancer binding factor 1 (ILF1)

Application Notes

Optimal dilution of the FOXK2 antibody should be determined by the researcher.

Immunogen

Amino acids TPPAAVREKGVQN were used as the immunogen for this FOXK2 antibody. This sequence is specific to isoform 1.

Storage

Aliquot and store the FOXK2 antibody at -20oC.