


FOXO1 Antibody (R35223)

Catalog No.	Formulation	Size
R35223-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2308
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This FOXO1 antibody is available for research use only.



Western blot testing of human umbilical cord lysate with FOXO1 antibody at 0.1ug/ml.
Predicted molecular weight ~70kDa.

Description

Additional name(s) for this target protein: Forkhead box O1

Application Notes

Optimal dilution of the FOXO1 antibody should be determined by the researcher.

Immunogen

Amino acids QDNKSSYGGMSQYNC were used as the immunogen for this FOXO1 antibody.

Storage

Aliquot and store the FOXO1 antibody at -20oC.