


FOXO3 Antibody (F42879)


Catalog No.	Formulation	Size
F42879-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F42879-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	O43524
Applications	Western Blot : 1:1000 IHC (Paraffin) : 1:10-1:50
Limitations	This FOXO3 antibody is available for research use only.



FOXO3 antibody immunohistochemistry analysis in formalin fixed and paraffin embedded human breast carcinoma.



FOXO3 antibody western blot analysis in HeLa, MCF-7 lysate. Expected molecular weight: 71-90 kDa depending on glycosylation level.

Description

This gene belongs to the forkhead family of transcription factors which are characterized by a distinct forkhead domain. This gene likely functions as a trigger for apoptosis through expression of genes necessary for cell death. Translocation of this gene with the MLL gene is associated with secondary acute leukemia. Alternatively spliced transcript variants encoding the same protein have been observed.

Application Notes

Titration of the FOXO3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 1-30 from the human protein was used as the immunogen for this FOXO3 antibody.

Storage

Aliquot the FOXO3 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.