


FOXO3A Antibody (R34068)

Catalog No.	Formulation	Size
R34068-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2309
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This FOXO3A antibody is available for research use only.



Western blot testing of human brain lysate with FOXO3A antibody at 0.1ug/ml. Expected molecular weight: 71-90 kDa depending on glycosylation level.

Description

Additional name(s) for this target protein: FKHRL1, AF6q21, FKHRL1P2; FOXO3

Application Notes

Optimal dilution of the FOXO3A antibody should be determined by the researcher.

Immunogen

Amino acids GAKQASSQSWVPG were used as the immunogen for this FOXO3A antibody.

Storage

Aliquot and store the FOXO3A antibody at -20oC.