



# FOXO4 Antibody (R33694)

Catalog No.	Formulation	Size
R33694-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

## **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	4303
Applications	Western Blot : 1-3ug/ml Enzyme ImmunoAssay : used as capture for recombinant protein (coat at 0.5-1ug/ml) ELISA (peptide) LOD : 1:64000
Limitations	This FOXO4 antibody is available for research use only.





Sandwich ELISA using FOXO4 antibody at 0.5ug/ml as the detect.

## **Description**

Additional name(s) for this target protein: MLLT7; myeloid/lymphoid or mixed-lineage leukemia (trithorax homolog, Drosophila) translocated to 7

## **Application Notes**

Optimal dilution of the FOXO4 antibody should be determined by the researcher.

## **Immunogen**

Amino acids EAASQDRMPQDLD were used as the immunogen for this FOXO4 antibody.

#### **Storage**

Aliquot and store the FOXO4 antibody at -20oC.