


## FOXP2 Antibody (R33432)


Catalog No.	Formulation	Size
R33432-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human, Rat
<b>Predicted Reactivity</b>	Cat, Cow, Dog, Mouse, Pig, Zebrafish
<b>Format</b>	Antigen affinity purified
<b>Host</b>	Goat
<b>Clonality</b>	Polyclonal (goat origin)
<b>Isotype</b>	Goat Ig
<b>Purity</b>	Antigen affinity
<b>Gene ID</b>	93986
<b>Applications</b>	Western Blot : 0.3-1ug/ml Immunohistochemistry (FFPE) : 2-4ug/ml ELISA (peptide) LOD : 1:64000
<b>Limitations</b>	This FOXP2 antibody is available for research use only.



Western blot testing of human cerebellum lysate with FOXP2 antibody at 0.3ug/ml. Predicted molecular weight: 80~83 kDa. Both bands block with immunizing peptide.



IHC staining of FFPE human prostate with FOXP2 antibody at 2.5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

## Description

Additional name(s) for this target protein: Forkhead box P2

## Application Notes

Optimal dilution of the FOXP2 antibody should be determined by the researcher.

1. This FOXP2 antibody will detect isoforms I, II, IV and V.

## Immunogen

Amino acids DEVEYQKRRSQKIT were used as the immunogen for this FOXP2 antibody.

## Storage

Aliquot and store the FOXP2 antibody at -20oC.