


FOXP2 Antibody (R33432)


Catalog No.	Formulation	Size
R33432-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Rat
Predicted Reactivity	Cat, Cow, Dog, Mouse, Pig, Zebrafish
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	93986
Applications	Western Blot : 0.3-1ug/ml Immunohistochemistry (FFPE) : 2-4ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This FOXP2 antibody is available for research use only.



Western blot testing of human cerebellum lysate with FOXP2 antibody at 0.3ug/ml.
Predicted molecular weight: 80~83 kDa. Both bands block with immunizing peptide.



IHC staining of FFPE human prostate with FOXP2 antibody at 2.5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Forkhead box P2

Application Notes

Optimal dilution of the FOXP2 antibody should be determined by the researcher.

1. This FOXP2 antibody will detect isoforms I, II, IV and V.

Immunogen

Amino acids DEVEYQKRRSQKIT were used as the immunogen for this FOXP2 antibody.

Storage

Aliquot and store the FOXP2 antibody at -20°C.