


## Foxp3 Antibody / Scurfin (RQ6868)

Catalog No.	Formulation	Size
RQ6868	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

**Bulk quote request**

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Mouse, Rat
<b>Format</b>	Antigen affinity purified
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit IgG
<b>Purity</b>	Antigen affinity purified
<b>Buffer</b>	Lyophilized from 1X PBS with 2% Trehalose
<b>UniProt</b>	Q99JB6
<b>Applications</b>	Western Blot : 1-2ug/ml
<b>Limitations</b>	This Foxp3 antibody is available for research use only.



Western blot testing of 1) rat thymus, 2) mouse spleen and 3) mouse thymus tissue lysate with Foxp3 antibody. Predicted molecular weight ~47 kDa.

## Description

FOXP3 (forkhead box P3) is a protein involved in immune system responses. The human FOXP3 genes contain 11 coding exons. Exon-intron boundaries are identical across the coding regions of the mouse and human genes. By genomic sequence analysis, the FOXP3 gene maps to the p arm of the X chromosome (specifically, Xp11.23). A member of the FOX protein family, FOXP3 appears to function as a master regulator in the development and function of regulatory T cells. While the precise control mechanism has not yet been established, FOX proteins belong to the forkhead/winged-helix family of transcriptional regulators and are presumed to exert control via similar DNA binding interactions during transcription.

## Application Notes

Optimal dilution of the Foxp3 antibody should be determined by the researcher.

## Immunogen

Amino acids EPEEFLKHCQADHLLDEK from the mouse protein were used as the immunogen for the Foxp3 antibody.

## Storage

After reconstitution, the Foxp3 antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.