


FRAT2 Antibody (R34119)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R34119-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|---------------------------|---|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Format | Antigen affinity purified |
| Host | Goat |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 23401 |
| Applications | Western Blot : 2-6ug/ml ELISA (peptide) LOD : 1:2000 |
| Limitations | This FRAT2 antibody is available for research use only. |



Western blot testing of human heart lysate with FRAT2 antibody at 2ug/ml. Predicted molecular weight ~25kDa.

Description

Additional name(s) for this target protein: Frequently rearranged in advanced T-cell lymphomas 2; FRAT2

Application Notes

Optimal dilution of the FRAT2 antibody should be determined by the researcher.

Immunogen

Amino acids GPDRIALQPSGSLL were used as the immunogen for this FRAT2 antibody.

Storage

Aliquot and store the FRAT2 antibody at -20oC.