


FRAT2 Antibody (R34119)

Catalog No.	Formulation	Size
R34119-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23401
Applications	Western Blot : 2-6ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This FRAT2 antibody is available for research use only.



Description

Additional name(s) for this target protein: Frequently rearranged in advanced T-cell lymphomas 2; FRAT2

Application Notes

Optimal dilution of the FRAT2 antibody should be determined by the researcher.

Immunogen

Amino acids GPDRIALQPSGSLL were used as the immunogen for this FRAT2 antibody.

Storage Aliquot and store the FRAT2 antibody at -20oC.				