


FRG1 Antibody (R36315)

Catalog No.	Formulation	Size
R36315-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2483
Applications	Western Blot : 2-4ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This FRG1 antibody is available for research use only.



Western blot testing of Jurkat nuclear lysate with FRG1 antibody at 2ug/ml.
Predicted/observed molecular weight: ~29kDa/35kDa.

Description

Additional name(s) for this target protein: Facioscapulohumeral muscular dystrophy region gene-1

Application Notes

Optimal dilution of the FRG1 antibody should be determined by the researcher.

Immunogen

Amino acids TVTNFGEISGT were used as the immunogen for this FRG1 antibody.

Storage

Aliquot and store the FRG1 antibody at -20oC.