


Frizzled 4 Antibody (R33413)

Catalog No.	Formulation	Size
R33413-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8322
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Frizzled 4 antibody is available for research use only.



Western blot testing of small intestine (Ileum) cell lysate with Frizzled 4 antibody at 1ug/ml

Description

Additional name(s) for this target protein: FZD4

Application Notes

Optimal dilution of the Frizzled 4 antibody should be determined by the researcher.

Immunogen

Amino acids DEEERRCDPIRIS were used as the immunogen for this Frizzled 4 antibody.

Storage

Aliquot and store the Frizzled 4 antibody at -20oC.