


Frizzled 8 Antibody (R36118)


Catalog No.	Formulation	Size
R36118-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8325
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : 3-5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Frizzled 8 antibody is available for research use only.



Western blot testing of human pancreas lysate with Frizzled 8 antibody at 1ug/ml.
Predicted molecular weight: ~73kDa.



IHC staining of FFPE human pancreas with Frizzled 8 antibody at 3ug/ml. HIER: microwaved with pH9 Tris/EDTA buffer, HRP-detect.

Description

Additional name(s) for this target protein: FZD8, FZ-8

Application Notes

Optimal dilution of the Frizzled 8 antibody should be determined by the researcher.

Immunogen

Amino acids SYPKQMPLSQV were used as the immunogen for this Frizzled 8 antibody.

Storage

Aliquot and store the Frizzled 8 antibody at -20oC.