


FRS2 Antibody (R36036)

Catalog No.	Formulation	Size
R36036-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10818
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This FRS2 antibody is available for research use only.



Western blot testing of 3T3-L1 lysate with FRS2 antibody at 0.1ug/ml. Predicted molecular weight: ~57kDa.

Description

Additional name(s) for this target protein: Fibroblast growth factor receptor substrate 2

Application Notes

Optimal dilution of the FRS2 antibody should be determined by the researcher.

Immunogen

Amino acids SRDEDDNLGPKTPS were used as the immunogen for this FRS2 antibody.

Storage

Aliquot and store the FRS2 antibody at -20oC.