


## FSD1 Antibody (R34684)

Catalog No.	Formulation	Size
R34684-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	79187
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This FSD1 antibody is available for research use only.



Western blot testing of human frontal cortex lysate with FSD1 antibody at 2ug/ml.  
Predicted molecular weight: ~55kDa.

## Description

Additional name(s) for this target protein: Fibronectin type III and SPRY domain containing 1

## Application Notes

Optimal dilution of the FSD1 antibody should be determined by the researcher.

## Immunogen

Amino acids KVQDIKAREKDGKGR were used as the immunogen for this FSD1 antibody.

## **Storage**

Aliquot and store the FSD1 antibody at -20oC.