


FSD1 Antibody (R34684)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R34684-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|-----------------------------|--|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Predicted Reactivity | Cow, Dog, Mouse, Rat |
| Format | Antigen affinity purified |
| Host | Goat |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 79187 |
| Applications | Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000 |
| Limitations | This FSD1 antibody is available for research use only. |



Western blot testing of human frontal cortex lysate with FSD1 antibody at 2ug/ml.
Predicted molecular weight: ~55kDa.

Description

Additional name(s) for this target protein: Fibronectin type III and SPRY domain containing 1

Application Notes

Optimal dilution of the FSD1 antibody should be determined by the researcher.

Immunogen

Amino acids KVQDIKAREKDGKGR were used as the immunogen for this FSD1 antibody.

Storage

Aliquot and store the FSD1 antibody at -20°C.